

Specifications

Applications	Cardiology/High Performance
Binaural Construction	Double Lumen
Brand	3M™ Littmann® Stethoscopes
Chestpiece Finish	Machined stainless steel
Chestpiece Technology	Single sided
Chestpiece Weight (metric)	90 g
Colour	Plum
Colour Group	Purple
Diaphragm Diameter (metric)	5.1 cm
Diaphragm Material	Epoxy/Fibreglass
Diaphragm Type	Tunable Diaphragm
Eartip Type	Soft Sealing
Environmental claim being made	Yes
Environmental claim relates to	Product
Evidence for Claim	Self Certified
Extra Eartips	No
Headset Material	Wide diameter aerospace alloy / Anodized aluminum
Model	Master Cardiology™
Occupation	Anesthesiologist, Cardiologist, Emergency Physician, EMT/EMS, Family Practitioner, Internist, Medical Student, Nurse, Pediatrician, Physician, Respiratory Specialist
Overall Length (Metric)	69 cm
Patient	Adult, Infant, Pediatric
Performance	10
Special Adaptors	Yes
Tube Colour	Plum