

# LABOUR DELIVERY RECOVERY BED LDR 100

EASY • ACCOMMODATING • SUPPORTIVE



## Technical Specification

Maximum height of the bed - 1150 mm

Minimum height of the bed - 610 mm

Overall Length 2130 mm

Overall Width with Siderails - 1020 mm

Mattress Size 900mm x 1950 mm x 125 mm

Maximum Head Angle from Horizontal 70°

Maximum back Angle (In Trendelenburg - Like Position) 7°

## Length of the mattress

Head section - 750 mm

Seat section - 545 mm

Leg Rest section - 655 mm

Thickness of the mattress - 125 mm

Sliding leg rest for comfortable care

## Standard Accessories

Catch basin, Mattress, Calf support

Hand grip, IV stand, Individual castors with brake  
hand set, Head board

## Optional Accessories

Full Leg Support



Call us toll free :

1 800 425 3503



## Phoenix Medical Systems Private Limited

No. 32/4, Jawaharlal Nehru Salai, Ekkattuthangal, Chennai 600 032, India.

Tel.: 91-44-2225 1177, 2225 1178 Fax : 91-44-2225 0194

✉ : sales@pmsind.com

🌐 : phoenixmedicalsystems.com