

Rotor Cat. No.	Description	
CR-17-24	Fixed angle rotor of CappRondo refrigerated centrifuge Capacity: 24 x 1.5mL/ 2.0mL microtubes 45°, Hole diameter 11.1mm, Max. height for tube fit 52mm Supplied with O-ring sealed	
CR-17-30	Fixed angle rotor of CappRondo refrigerated centrifuge Capacity: 30 x 1.5mL 2.0mL microtubes 45°, Hole diameter 11.1mm, Max. height for tube fit 52mm Supplied with O-ring sealed	
CR-17-18	Fixed angle rotor of CappRondo refrigerated centrifuge Capacity: 18 x 1.5mL/ 2.0mL microfilter tubes 45°, Hole diameter 11.1mm, Max. height for tube fit 65mm Supplied with O-ring sealed	
CR-17-5	Fixed angle rotor of CappRondo refrigerated centrifuge Capacity: 12 x 5mL screw and snap cap microtubes 45°, Hole diameter 17mm, Max. height for tube fit 64mm Supplied with O-ring sealed	

<p>CR-17-PCR</p>	<p>Fixed angle PCR-tube rotor of CappRondo refrigerated centrifuge Capacity: 64 x 0.2mL PCR tubes or 8 rows of PCR strips 45°, Hole diameter 6.5mm, Max. height for tube fit 16mm</p>	
<p>CR-17-D8</p>	<p>Drum Rotor of CappRondo refrigerated centrifuge Capacity: 8 x 1.5mL/ 2.0mL microtubes (buckets are not included) 90°, Hole diameter 11.1mm, Max. height for tube fit 50mm</p>	
<p>CR-17-D16</p>	<p>Drum Rotor of CappRondo refrigerated centrifuge Capacity: 16 x 1.5mL/ 2.0mL microtubes (buckets are not included) 90°, Hole diameter 11.1mm, Max. height for tube fit 43mm</p>	