

POINTCARE M3 REAGENT DISC DETAIL

DISC	TEST
General Chemistry I	TP ALB GLO* A/G* TBIL DBIL IBIL* ALT AST TG CHOL HDL-C LDL-C GLU CRE UREA UA
Clinical Emergency	AST CK CK-MB LDH α-HBDH GLU AMY CRE UA tCO ₂ K ⁺ Na ⁺ Cl ⁻
Renal Function	ALB CRE UREA UA tCO ₂ Ca ²⁺ PHOS
Liver Function	TP ALB GLO* A/G* TBIL DBIL IBIL* ALT AST GGT ALP
Myocardial Enzyme	AST CK CK-MB LDH α-HBDH
Electrolyte	tCO ₂ Ca ²⁺ PHOS K ⁺ Na ⁺ Mg ²⁺ Cl ⁻
Glucose & Lipid	TG CHOL HDL-C LDL-C GLU GSP
GLU, Lipid & HCY#	TG CHOL HDL-C LDL-C GLU HCY#
General Chemistry II	GLU AMY CRE UREA tCO ₂ K ⁺ Na ⁺ Cl ⁻
Liver & Renal	TP ALB GLO* A/G* TBIL ALT AST GGT GLU CRE UREA
CRP & SAA#	CRP# SAA#