

DISPOSABLE SAFETY FACE SHIELDS

Fluid Resistant Full-Face Visor

For the protection of the facial area and associated mucous membranes (eyes, nose, mouth) from splashes, sprays, and spatter of body fluids.

Key Specifications:

- Optically clear
- No distortion
- Wrap around
- Fluid Resistant
- Full Face Visor

Ideal for:

- Hospital/Primary Care
- Allied Health Professionals
- Aged Care Workers
- Workplace H&S
- Retail workers

Face shields provide a barrier to acutely-expelled aerosols of body fluids and are commonly used as an alternative to goggles as they confer protection to a larger area of the face.



DISPOSABLE SAFETY FACE SHIELDS

Fluid Resistant Full-Face Visor

Key Specifications cont:

Visor/Shield:

- Sized to: 33cm x 22cm
- 2cm elastic Band - domed to each side
- Adhesive Tape on rear to apply Foam
- PET of 0.285 anti fog I Clear
- Supplied with protective film both sides "to be removed

Foam Band:

- Sized to: 23cm x 3.5cm @ 3.3cm depth
- Density: 8-20g / cm²
- Colours - White | Blue | Black



Presented: 10pcs Foam | 10pcs Visor - in separate sealed poly bag



**Due to Covid 19 Demand/
Stock shortages - Only Black
available until April 15th**

