

PRODUCT CODE:
65976

PRODUCT NAME:
Help-It General Purpose Ear Loop Mask Child Size
(Box = 50 Masks / Carton = 40 Boxes)



DESCRIPTION

Provides filtration to reduce the amount of inhaled particulate matter. Fluid resistant and provides the wearer protection against large droplets, splashes, or sprays of bodily or other hazardous fluids. Provides source control by containing the wearer’s respiratory secretions, droplets, and aerosols.

The Help-It® face mask for consumers covers the mouth, nose and chin and is made of non-woven fabric. The elastic ear loops allow the face mask to be easily put on and removed and are sufficiently robust to hold the face mask in place. Single strand plastic coated nose wire allows the face mask to be molded into place around the nose.

SPECIFICATIONS

Small size 145mm x 95mm
EN14683:2019+AC:2019 ANNEX 3
EN14683:2019 Differential Pressure <6mmH20/cm²
BFE ≥ 98%
ISO22609:2004 Splash Resistance 160mmHg

FEATURES

Disposable, light weight, fluid resistant, latex free, comfortable, breathable, manufactured in China.

APPLICATIONS

Public transport, air travel, shopping, schools, events, hospitality, public gatherings.

RESPONSIBLE DISPOSAL

Do not handle the front of the face mask. Discard in appropriate secure rubbish bin. Wash hands with soap and water or rub them with hand sanitizer.

STORAGE

Store in a cool, dry place out of direct sunlight.

EXPIRATION DATE

5 years shelf life from date of manufacture.

LIMITATIONS

Face masks should not be placed on very young children (<2 years) or anyone who has trouble breathing or unable to remove the face mask without assistance. The face mask contains metal and should not be worn during MRI procedures.